

# Join The Community RADRoller.com

#### Lifetime Guarantee.

Lifetime guarantee against manufacturer defects, barring misuse, prolonged sun exposure, anticipated wear, and loss. Contact **Service@RADRoller.com** for more information.

### Legal:

Use at your own risk. **RAD Helix** is not meant to treat, cure, or prevent any disease or condition. Talk to your doctor before using the **RAD Helix**. Please use responsibly.

Patent Pending.











Once you get your hands on a RAD Helix, you're going to want to get started right away. Hey, we get it. These things are, well, rad. But before you dive into the following (and your own) RAD movements, there are a few key things you need to know.\*

- Only keep pressure on a given location for 7–10 seconds.
  Any longer and you can cut off the muscle's blood supply, which can result in soreness. No good.
- Roll both lengthwise and across the muscle. The power of the **RAD Helix** is in its ability to break through restrictions in three dimensions, so take advantage.
- Don't roll directly on bone. That's just plain irresponsible.



The bones of the low back bear the entire weight of the spine, so let's make sure they're aligned and happy. Humans are meant to have a natural curve in the low back that can slowly deteriorate from prolonged sitting or flexion (like cycling). The RAD Helix is meant to open up these curves to provide extension, while giving you the ability to flush out toxins and break up tough tissue in the entire spine. Used in conjunction with the RAD Roller, you'll never be without a healthy release for your back.





#### **Low Back**

To help open up the lower back, simply lay back with the roller directly beneath your belly button. Relax and breathe naturally. Move up the back towards your head, stopping and relaxing at any tight spots.

## **Neck and Upper Back**

The **RAD** Helix can provide a gentle release for your neck, and a primer for using the **RAD** Roller or **RAD** Rounds. Simply lay back with the roller perpendicular to your neck. Raise your knees up, and gently push back and forth with your feet to provide a gentle massage to the superficial muscles of the neck and upper back. The **RAD** Helix is also a perfect primer or warm up for digging deeper with the **RAD** Roller.